C10RN

Convex Black&White

www.konted.cn

KONTED



Frequency	Depth	Mode	Applications
3.5-5MHz	90-200mm	B Mode B+M	Abdominal Emergency Cardiac Urology Obstetrics / Gynecology, Lung Kidney

Download our ultrasound software from App Store or Google Play



- \checkmark Unlimited users
- ✓ Free software
- ✓ No hidden fees
- ✓ Support printing reports
- ✓ Chinese, English, Russian, ItalianSpanish, Portuguese (Brazil)

Software Controls GN(gain) D(depth) ENH(enhancement) DR(dynamic range) F(frequency) Mode 8TGC Biopsy Annote Page left/right up/down Google play

Measurement


B: Length, Area/Circumference,

Angle,Trace,Distance,GA(CRL,BPD,GS,FL,HC,AC) EFW(BPD,FL)

B+M: Heart Rate, Time, Distance

Offer great images to everyone





Package includes:

- 1* wireless probe
- 1* USB charger
- 1* plug

- 1* wireless charger
- 1* wrist strap
 - 1* portable aluminium alloy box

Connection Type

WiFi (802.11n/20MHz/5G) Probe built-in WiFi signal

Supporting Devices Tablet, Smart phone, PC, etc

Working system Android,IOS,Windows

Buttons

3 programmable buttons Lithium Battery $4200 \text{ mAh. } 3.7 \vee$ $4000 \text{ mAh. } 3.7 \vee$ $4000 \text{ mAh$

Net weight235gElements128FootprintR60Boot-up



Platform dependent, generally <5 seconds



- * Training how to install the machine, training how to use the machine.
- * Spare parts replacement and other software/hardware trouble shooting
- * Offer OEM&ODM service for distributors
- * Offer SDK depends on different models
- * Free software and free upgrade the software
- * 18 Months Warranty, extendable up to 3 years

BEIJING KONTED MEDICAL TECHNOLOGY CO.,LTD Tel:8610-60219113 Fax:8610-60219213 Service and Cooperation :info@konted.cn Website:www.konted.cn Facebook: Beijing Konted Youtube: Medical Konted Address: Room 111,1F, Building 3, No. 27 Yongwang Road,Daxing Biological Pharmaceutical Industry Base,Daxing District,Beijing 102629,China.