



2 IN 1 DUAL HEAD

192 Elements 64 Channels HD ULTRASOUND

Ultrasound

C10XS CONVEX & LINEAR PROBE



KONTED

www.kontedmed.com

Probe Specs

Model: C10XS
Convex Freq: 3.2/5Mhz
Scan Depth: 90-305mm
Linear Freq: 7.5/10Mhz
Scan Depth: 20-100mm
Elements: 192
Channel: 64
Waterproof: IP68
Mode: B.B+B.B+M
Color.PW. PDI.

Size

Dimensions:
160x65x18mm
Weight:263g
Radius:R60 (Convex)
Footprints:L45 (Linear)

14+ Presets

Abdomen
GYN&OB
Cardiac, Urology,
Kidney, Lung

Thyroid
Small parts
pediatrics, Vascular
Carotid, Breast, Msk
Nerve

Image Management

JPEG, MP4, Dicom
Save&Share
image/video/report


Connection (WiFi)

WiFi:
802.11n/20MHz/2.4G/5G
System:
Android. IOS. Windows

Device:
Android :Tablet/Phone/TV
IOS :Ipad/iPhone (only WiFi)
Windows PC

Battery

Capacity: 2800 mAh
Working time: 3h
48h Ultra-long standby
Recharge time: 2~3h by
Wireless charging



Abdominal & Superficial

FDA CE
ISO TGA